

Meet Our Extended Family


SKATER[™]
COMPLETE
DRAINAGE SYSTEM

The SKATER[™] All Purpose and Nephrostomy Drainage Catheter portfolio includes longer and shorter lengths to meet the needs of any patient.


ARGON
MEDICAL DEVICES

SKATER™

COMPLETE DRAINAGE SYSTEM



SKATER All Purpose and Nephrostomy Locking (Box of 5)

CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET	TROCAR SIZE
756606015	6F x 15cm	.035"	Yes	19G
756606020	6F x 20cm	.035"	Yes	19G
756606025	6F x 25cm	.035"	Yes	19G
756606035	6F x 35cm	.035"	No	
756607015	7F x 15cm	.035"	Yes	18G
756607020	7F x 20cm	.035"	Yes	18G
756607025	7F x 25cm	.035"	Yes	18G
756608015	8F x 15cm	.038"	Yes	17G
756608020	8F x 20cm	.038"	Yes	17G
756608025	8F x 25cm	.038"	Yes	17G
756608030	8F x 30cm	.038"	Yes	17G
756608035	8F x 35cm	.038"	No	
756608045	8F x 45cm	.038"	No	
756610015	10F x 15cm	.038"	Yes	17G
756610020	10F x 20cm	.038"	Yes	17G
756610025	10F x 25cm	.038"	Yes	17G
756610030	10F x 30cm	.038"	Yes	17G
756610035	10F x 35cm	.038"	No	
756610045	10F x 45cm	.038"	No	
756610060	10F x 60cm	.038"	No	
756612015	12F x 15cm	.038"	Yes	14G
756612020	12F x 20cm	.038"	Yes	14G
756612025	12F x 25cm	.038"	Yes	14G
756612030	12F x 30cm	.038"	Yes	14G
756612035	12F x 35cm	.038"	No	
756612045	12F x 45cm	.038"	No	
756612060	12F x 60cm	.038"	No	
756614015	14F x 15cm	.038"	Yes	14G
756614020	14F x 20cm	.038"	Yes	14G
756614025	14F x 25cm	.038"	Yes	14G
756614030	14F x 30cm	.038"	Yes	14G
756614035	14F x 35cm	.038"	No	
756614045	14F x 45cm	.038"	No	
756614060	14F x 60cm	.038"	No	
756616025	16F x 25cm	.038"	Yes	14G
756616030	16F x 30cm	.038"	Yes	14G
756616035	16F x 35cm	.038"	No	
756616045	16F x 45cm	.038"	No	
756616060	16F x 60cm	.038"	No	

*Flexible stiffening cannula not available on 6F and 7F.



© 2021 Argon Medical Devices, Inc. All Rights Reserved.

No reproduction of any portion of this work shall be made without prior permission of the author.

SKATER All Purpose and Nephrostomy Non-Locking (Box of 5)

CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET	TROCAR SIZE
756106015	6F x 15cm	.035"	Yes	19G
756106020	6F x 20cm	.035"	Yes	19G
756106025	6F x 25cm	.035"	Yes	19G
756107020	7F x 20cm	.035"	Yes	18G
756107025	7F x 25cm	.035"	Yes	18G
756108020	8F x 20cm	.038"	Yes	17G
756108025	8F x 25cm	.038"	Yes	17G
756108030	8F x 30cm	.038"	Yes	17G
756110020	10F x 20cm	.038"	Yes	17G
756110025	10F x 25cm	.038"	Yes	17G
756110030	10F x 30cm	.038"	Yes	17G
756112020	12F x 20cm	.038"	Yes	14G
756112025	12F x 25cm	.038"	Yes	14G
756114020	14F x 20cm	.038"	Yes	14G
756114025	14F x 25cm	.038"	Yes	14G
756116025	16F x 25cm	.038"	Yes	14G

*Flexible stiffening cannula not available on 6F and 7F.

SKATER Biliary Drainage Catheter w/ Locking Pigtail (Box of 5)

CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE
755308040	8F x 40cm	.038"
755310040	10F x 40cm	.038"
755312040	12F x 40cm	.038"

SKATER Biliary Drainage Catheter w/ Non Locking Pigtail (Box of 5)

CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE
755208040	8F x 40cm	.038"
755210040	10F x 40cm	.038"
755212040	12F x 40cm	.038"

Each Set Includes

- (1) SKATER Catheter
- (1) Removable pigtail straightener
- (1) Metal Stiffening Cannula
- (1) Flexible Stiffening Cannula (not available on 6F and 7F)
- (1) Locking system with snap-off tab (not available on non-locking pigtail catheters)



96-6400-53/0521A